

Hybrid
White Knight Implants LLC




Patient presents at initial visit.



Transfers connected with duralay in preparation for master impression.



Master Impression with analogs in place on transfers.




Temporary copings in place on master soft tissue model.



Verification Jig on master model.
Temporary cylinders connected
with triad gel.



Esthetic setup ready for try in.




Master model, setup and verification jig are sent to 3I for bar fabrication. Virtual images are e-mailed to the lab for final design ok.



Finished titanium bar is returned
to the laboratory.







Surface is sandblasted to increase retention.



Pink opaque is applied to mask bar
beneath acrylic.



Wax up tried in on bar for final
check of fit and esthetics.




Processed Hybrid denture.



Screw holes covered with cotton pellets
and composite.





Book design by Howard Gershman
White Knight Implants LLC