



**Cam StructSURE Overdenture
White Knight Implants LLC**





Open tray and soft tissue master model.



Anchored bite block and verification stent.




Verification stent tried in the patient's mouth.



Esthetic setup.



This step is tried in the mouth and sent to 3i along with the master model and the verification stent.




At the factory technicians scan the setup and send a virtual image of the bar to the laboratory for approval.



The packaged bar is returned to the laboratory.



Note the accurate fit of the bar on the master model.



Retentive elements are placed on the bar to be incorporated into the setup.





The bar is tried in the mouth.






The reconstituted setup in place on the bar.



Before processing a mesh framework is waxed on a refractory model, cast and opaqued with pink light -cured material.







The finished case on the master cast.





Clips and opaqued mesh shown with the finished prosthesis.



Finished case seated in the mouth.



Restorative and Intra-oral Photography: Gary Graham DMD
Prosthetics: White Knight Implants LLC

11 Research Drive 3B
Woodbridge, CT 06525
203-553-9990
www.WhiteKnightImplants.com